



DO106024

## Lateral sliding holder

Accessory to move a component from the optical axis.

### School level

Senior high school.

### Technical features

Slab put up on sliding holder for prismatic bench  
Adjustment range: 50 mm (+/- 25mm)  
Graduation: 1.0 mm  
Adjustment: 1mm/tour  
Linear guidance with pads  
Dimensions: 130 x 100 x 50 mm  
Weight: 250 g.  
Packaging: individual box

### Product advantages

- Adjustment with endless crew for **high precision**.
- **Easy handling** thanks to the knurled button.

### Associated products

Prismatic premium optics bench – **DO106023**  
Diffraction suitcase – **DO115021**  
Laser TRIO red (also available in green) – **DO108015**  
Slide 4 diffraction holes – **DO115015**  
Slide 6 simple slits – **DO115018**  
Slide 7 slits 7 lines – **DO115022**  
Slide double slits – **DO115014**  
Slide triple slits of Young – **DO115009**



*Utilisation sur BO prismatique*